


Scanning electron microscopy & microtomic techniques

2017.1.30


Lab meeting

Yu-Ling Huang

Resolving power of microscopes



	Magnification	Resolution	Depth of field
Light microscope	1-1000x	200nm	2μm
SEM	2-300,000x	3.5-6nm	Up to 1 cm
TEM	2,00-900,000x	0.21-0.5nm	0.11μm



$$r = \frac{1.22\lambda}{2n \sin \theta} = \frac{0.61\lambda}{NA}$$

where

r is the minimum distance between resolvable points, in the same units as λ is specified

λ is the wavelength of light, emission wavelength, in the case of fluorescence,

n is the index of refraction of the media surrounding the radiating points,


θ is the half angle of the pencil of light that enters the objective, and

NA is the numerical aperture

Visible light wavelength= 400-700nm

Electron wavelength= 0.004nm

Electron microscope



SEM

- Useful for observing surface details
- Electron beam is focused on surface of specimen with electromagnetic lenses in a vacuum chamber
- Electron beam scans back and forth across specimen
- Electron beam emit different electrons from the sample surface






FIG. 2.24 An illustration of how the scanning electron microscope can reveal surface relief when used with a secondary electron detector

SEM sample preparation

- Fixation
- Dehydration
- Critical point drying
- Coating



Fixation

- Problems of biological specimens:
 - Contrast lost
 - Redistribution or degradation of organelles
 - Wound induced changes
- Fixation preserves sample morphology, anatomy, cell and subcellular structure to keep the sample close to living status
- **Chemical fixation:** slow, general method.
- **Cryofixation:** fast, but ice crystal may break the tissue. Shrinkage is another problem.
- **Microwave fixation:** fast, but heat may affect subcellular structure.

Factors affect chemical fixation

- sample size and penetration rate
- pH and osmolarity
- Freshness and concentration of fixatives
- Fixation duration
- Temperature
- Amount of fixative
- Air space in tissue



Fixatives

- Coagulating Fixatives:
 - fix rapidly, change hydration state, proteins are coagulated. Not for high resolution
 - Example: acetone, methanol, ethanol
- Non-coagulating fixative
 - Slow penetration, form crosslink to proteins. Good for high resolution
 - Example: formaldehyde, glutaraldehyde, OsO_4

General fixatives for SEM

- 2.5-4% Glutaraldehyde (GA) in phosphate buffer (pH 6.8-7.4)
 - GA is good for preserve ultra-structure
 - Slow penetration
- Karnovsky's fixative (half-strength)
 - 2% Paraformaldehyde (fast)
 - 2.5% Glutaraldehyde (slow)
 - 0.1M phosphate Buffer (pH 6.8-7.4)

Dehydration




- SEM samples are in vacuum chamber, so the samples must be 100% dry
- Dehydration in ethanol series to remove water from tissue through substitution without causing extreme shrinkage
- Dehydration time depends on sample size, but general is 10-20min.
- Do not store samples in ethanol for extended time which may cause shrinkage (except in 70% ethanol)
- For samples with high water content, starts from 30% ethanol, otherwise 50% ethanol is fine.
- Keep samples 'wet' all the time

(30%→)50%→70%→80%→90%→95%→100%

Critical point drying (CPD)



- Water leaves tissue could cause breakage due to high surface tension
- No surface tension, no liquid/gas boundary in supercritical fluid
- CO₂ is safe and has lower temperature and pressure than water to reach critical point
- Ethanol or acetone is compatible with CO₂



Critical point dryer



Sample mounting



Carbon adhesive tapes/ tabs



SEM sample stub

Sputter coating

- Biological samples is non-conductive
- Protect samples from heat damage
- Reduce sample charging
- Gold or platinum
- Thickness
- Direction



Hitachi S-4800 Type II Ultra-High Resolution Field Emission Scanning Electron Microscope




Microtomy


- Microtome: a tool for cutting samples into thin slices
- Normal thickness of the section: 1-10 μm
- Ultra-thin section: $< 1 \mu\text{m}$ (TEM)
- Types of microtome:
 - Free hand sections
 - Sliding/sledge microtome
 - Rotary microtome
 - Cryomicrotome
 - Ultramicrotome

Free hand section


- Fast
- Thick sections
- Temporary preparations





Free hand section







Casuarina equisetifolia L. (Casuarinaceae)






Sliding/sledge microtome



<https://youtu.be/VcOmtP566b4>

Rotary microtome



<https://youtu.be/qGyE4p9XanA>

Cryomicrotome



Ultramicrotome



Microtome knife

- Steel knife





Microtome knife

- Steel knife
- Razor blade




Microtome knife

- Steel knife
- Razor blade
- Disposable steel blade




Microtome knife

- Steel knife
- Razor blade
- Disposable steel knife
- Glass knife




Making Glass Knives




1. Clean prepared glass strips
2. Score at right angles with glass knife maker



3. Break into squares (25mm)



4. Score diagonally and break into two knives



Microtome knife

- Steel knife
- Razor blade
- Disposable steel blade
- Glass knife
- **Diamond knife**



Embedding medium



Paraffin

- Low cost
- Fast
- Soft and easy to cut
- Need de-paraffin before staining
- Section 5-10 μm




Technovit 7100 (Glycol Methacrylate)

- High cost
- Slow
- Hard to cut
- Easy to stain
- Section 1-5 μm

Sample preparation for rotary microtome

- Fixation
- Dehydration
- Infiltration
- Embedding



Source: Mescher AL: Junqueira's Basic Histology: Text and Atlas, 12th Edition: <http://www.accessmedicine.com>
Copyright © The McGraw-Hill Companies, Inc. All rights reserved.

Fixation

- Similar to the SEM process
- Common fixative in botanical studies:
 - Formalin- Acetic Acid –Alcohol (FAA)

50% (or 70%) Ethanol	90 ml
Glacial Acetic Acid	5 ml
Formaldehyde (38%)	5 ml

Dehydration

- Select the dehydration solution which is compatible with the embedding medium
- Series of t-butanol (TBA) in 95% ethanol and water for paraffin
- Series of Ethanol for Technovitt 7100
50%→70%→80%→90%→95%→100%




TBA dehydration series

steps	t-butanol	95% ethanol	dH ₂ O
1	10	40	50
2	20	50	30
3	35	50	15
4	55	45	0
5	75	25	0
6	100	0	0

Paraffin infiltration


- Let embedding medium infiltrate into tissue
- Paraffin infiltration in oven around 60 °C
- Lower temperature is better but must higher than the melting point of paraffin



Technovit 7100 infiltration

Prep	100% Ethanol	“Technovit 7100”
I	5	1
II	3	1
III	2	3
IV	1	5


“Technovit 7100”: 100mL Technovit 7100 + 1 g hardener I (1 pack)



Technovit 7100 infiltration

Prep	100% Ethanol	“Technovit 7100”
I	5	1
II	3	1
III	2	3
IV	1	5

“Technovit 7100”: 100mL Technovit 7100 + 1 g hardener I (1 pack)



Infiltration steps:

Technovit 7100 Prep I 12 hours

Technovit 7100 Prep II 12 hours

Technovit 7100 Prep III 12 hours

Technovit 7100 Prep IV 12 hours

“Technovit 7100” 12 hours

“Technovit 7100” 24 hours**

**If longer than 24 hours, put in 4 °C to prevent polymerization.




Paraffin embedding



Paraffin sample block

- Trimming
- Attach to a wood block



Technovit 7100 embedding

- Prepare embedding medium:


15 mL “Technovit 7100”

1 mL Technovit 7100 hardener II

0.6 mL PEG 400

- Fill each block of HistoForm S with the embedding medium using micropipette (P1000)
- Transfer each sample to each block with tweezers, toothpicks or pipette tips
- Wait samples sink to bottom, and adjust the orientation of the samples
- Put the HistoForm S horizontally in 4 °C, overnight (~ 1 day) for polymerization


Technovit 7100 sample block



*Put sample block on 80°C heat plate to soften the block before trimming


** re-embedding if the sample orientation is not right

Microtome operation



AO Spencer 815 Rotary Microtome

Knife angle (clearance angle)




paraffin section: 4-6 degree

Technovit : 10 degree


Paraffin ribbon



Paraffin sections could form a straight and long ribbon if everything is good



Inadequate
trimming

Technovit 7100 section





Staining



Turkish Journal of Botany

<http://journals.tubitak.gov.tr/botany/>

Research Note

Turk J Bot

(2013) 37: 784-787

© TUBİTAK


doi:10.3906/bot-1204-9

New staining technique for fungal-infected plant tissues

João Paulo Rodrigues MARQUES, Marli Kasue Misaki SOARES, Beatriz APPEZZATO-DA-GLORIA⁺
Department of Biological Sciences, Luiz de Queiroz College of Agriculture, University of São Paulo, Piracicaba, Brazil

2.4. Staining procedure

1. Staining in 5% blue cotton in lactophenol for 20 min (Macedo, 1997);
2. Three 1-min washes in distilled water;
3. Staining in 1% aqueous safranin O for 10 s;
4. Three 1-minute washes in distilled water;
5. After drying, the slide can be mounted with an Entellan[®] synthetic resin.



Useful references

- Sass, JE. 1958. Botanical Microtechnique. 3rd Ed. Iowa State University Press, Ames, Iowa.
- Richerds, OW. 1942. The effective use and proper care of the microtome. Spencer Lens Company, Buffalo, NY.
- Yeung, EC and Chan, CKW. 2015. **Plant Microtechniques and Protocols..** Springer International Publishing, Switzerland.
- Yeung, EC and Chan, CKW. 2015. Glycol methacrylate: the art of embedding and serial sectioning. Botany 93: 1–8 (2015)
dx.doi.org/10.1139/cjb-2014-0177
- Dykstra, MJ. 1993. A manual of applied techniques for biological electron microscopy. Plenum Press, New York, NY.
- Leica Biosystems Education page
<http://www.leicabiosystems.com/pathologyleaders/>